Jpn. J. Protozool., 39(1), 2006 pp. 113-114 Japan Society of Protozoology

Loop-mediated isothermal amplification (LAMP) and its application for the detection of African trypanosomes

Noboru INOUE, Oriel M. M. THEKISOE, Noritaka KUBOKI, Chihiro SUGIMOTO (Obihiro University of Agriculture and Veterinary Medicine, National Research Center for Protozoan Diseases, Research Unit for Advanced Preventive Medicine)

While PCR is a method of choice for the detection of African trypanosomes in both man and animals, the expense of this method negates its use as a diagnostic method for the detection of African trypanosomosis in African countries where it is endemic. The loop-mediated isothermal amplification (LAMP) reaction is a method that amplifies DNA with high specificity, efficiency and rapidity under isothermal conditions using only simple incubators. An added advantage of LAMP over PCR-based methods is that DNA amplification can be monitored spectrophotometrically and/or with the naked eye, without the use of dyes. Here we present our conditions for a highly sensitive, specific and easy diagnostic assay based on LAMP technology for the detection of parasites of the *Trypanosoma brucei* group (including *T. b. brucei*, *T. b. gambiense*, *T. b. rhodesiense*, and *T. evansi*) and *T. congolense*. We show that the sensitivity of the LAMP-based method in detecting trypanosomes in vitro is up to 100 times higher than that of PCR-based methods in laboratory conditions. *In vivo* studies in mice infected with human-infective *T. b. gambiense* further highlight the potential clinical importance of LAMP as a diagnostic tool for identification of African trypanosomiasis.